

How to use the EpiPen® / EpiPen® Jr. (Epinephrine) Auto-Injectors



Remove yellow or green cap from carrying case

- Grasp the Auto-Injector with black tip pointing downward
- Pull off grey safety cap



Place black tip against mid-outer thigh and press firmly until the Auto-Injector activates. Hold while counting for several seconds, then remove

- Massage the injected area for 10 seconds



After administration

Call 911 or have someone take you to the emergency room

After using EpiPen® / EpiPen® Jr. follow 3 easy safety steps:



Carefully place used Auto-Injector, needle-end first, into storage tube



Screw cap of carrying case on completely

- This automatically bends needle back and secures pen so it won't fall out of tube



Give any used Auto-Injectors to emergency responders or emergency room personnel



EpiPen Auto-Injectors are indicated for the emergency treatment of anaphylactic reactions and for patients determined by a physician to be at increased risk for anaphylaxis.

For additional information, please visit www.EpiPen.ca